
EN 3	<b>Biodesign® Anal Fistula Plug</b> Instructions for Use
CS 5	<b>Biodesign® Anální píštělová zátka</b> Návod k použití
DA 8	<b>Biodesign® anal fistelprop</b> Brugsanvisning
DE 11	<b>Biodesign® Analfistelverschluss</b> Gebrauchsanweisung
EL 14	<b>Biodesign® Βύσμα περιεδρικού συριγγίου</b> Οδηγίες χρήσης
ES 18	<b>Biodesign® Tapón para fístulas anales</b> Instrucciones de uso
ET 21	<b>Biodesign® Päraku-uurise kork</b> Kasutusjuhend
FR 24	<b>Biodesign® Obtuteur de fistule anale</b> Mode d'emploi
HU 27	<b>Biodesign® anális fisztuladugó</b> Használati utasítás
IT 30	<b>Biodesign® Tampone per fistola anale</b> Istruzioni per l'uso
NL 33	<b>Biodesign® anale fistelplug</b> Gebruiksaanwijzing
NO 36	<b>Biodesign® analfistelplugg</b> Bruksanvisning
PL 39	<b>Biodesign® Zatyczka przetoki odbytniczej</b> Instrukcja użycia
PT 42	<b>Biodesign® Plug para fístula anal</b> Instruções de utilização
SV 45	<b>Biodesign® analfistelplugg</b> Bruksanvisning

\*FP0037-03K\*

## BIODESIGN® ANAL FISTULA PLUG

### INTENDED USE

The Cook® Biodesign® Anal Fistula Plug is for implantation to reinforce soft tissue where a rolled configuration is required, for repair of anorectal fistulas. The plug is supplied sterile and is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**ANAL FISTULA PLUG** This symbol means the following: Anal Fistula Plug

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

- This plug is derived from a porcine source and should not be used in patients sensitive to porcine materials.
- **Not for vascular use.**

### PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- Plug is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard plug if mishandling has caused possible damage or contamination, or if the plug is past its expiration date.
- Do not implant the plug in a grossly infected or abscessed fistula tract. **A seton should be used until there is no evidence of acute inflammation, purulence, or excessive drainage. Allow the tract to mature and stabilize for six to eight weeks before placing the plug.**
- Ensure that the plug is rehydrated prior to placement, cutting, or suturing.

### GENERAL

- Users should be familiar with surgical technique for anorectal fistula repair.
- Users should exercise good surgical practice for the management of clean-contaminated, contaminated or infected fields.
- The potential for infection of the graft material following implantation may be reduced by cleaning of the fistula tract.
- The plug should be handled aseptically, minimizing contact with latex gloves.
- The plug should be introduced through the internal (primary) fistula opening.
- The plug should be drawn into the fistula tract only until the internal (primary) opening is satisfactorily closed / occluded.
- The external (secondary) opening of the fistula tract should remain open.
- **IMPORTANT:** Users should counsel patients on abstaining from strenuous activity for a period of six weeks after plug placement. See Post-Operative Care section.
- **Expect some drainage for two to four weeks. Drainage can continue for up to 12 weeks after the procedure as the plug is incorporated and the fistula tract is closed.**

### POTENTIAL COMPLICATIONS

Complications that can occur with the plug include:

- |   |                      |             |
|---|----------------------|-------------|
| • Inflammation                                | • Induration         | • Migration |
| • Extrusion                                   | • Seroma formation   | • Infection |
| • Abscess                                     | • Fistula recurrence | • Sepsis    |
| • Swelling                                    | • Redness            | • Bleeding  |
| • Delayed or failed incorporation of the plug |                      |             |

If any of the following conditions occur and cannot be resolved, plug removal should be considered:

- Infection
- Abscess
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

### STORAGE

This plug should be stored in a clean, dry location at room temperature.

### STERILIZATION

This plug has been sterilized with ethylene oxide.

### INSTRUCTIONS FOR USE

#### Required Materials

- A sterile dish (kidney dish or other bowl)

- Rehydration fluid: room temperature sterile saline or sterile lactated Ringer's solution.
- Suitable resorbable suture, such as: 0 chromic, 2-0 or 0 coated polyglycolic acid suture (coated PGA).
- 10cc syringe
- A short 14ga or 16ga catheter
- Saline or hydrogen peroxide for flushing

**NOTE:** Always handle the plug using aseptic technique, minimizing contact with latex gloves.

#### PREPARATORY

1. Remove the packaging containing the plug from the box.
2. Remove the inner pouch containing the plug from the outer package using aseptic technique. Place the inner pouch in the sterile field.
3. Using sterile gloved hands, open the inner pouch carefully, and aseptically remove the plug with a sterile instrument. Place the plug into the sterile dish in the sterile field.
4. Add enough rehydration fluid to the dish to fully submerge the plug. Allow the plug to rehydrate, fully submerged, until the desired handling characteristics are achieved, but do not allow the plug to rehydrate any longer than 2 minutes.
5. Attach a suitable resorbable suture (approximately 30 cm in length) around the tail (narrow end) of the plug for pulling it into the fistula tract.
6. Prepare the patient and surgical site using standard surgical techniques appropriate for anal fistula repair.

**NOTE:** The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing through the use of enemas and cathartic agents, and administration of prophylactic antimicrobial agents, oral or intravenous. Insufficient cleansing or inadequate antibacterial prophylaxis can predispose the patient to infections.<sup>1-3</sup>

#### PROCEDURAL


1. Perform under local, regional, or general anesthesia.
2. If a seton is not already in place, identify the internal (primary) fistula opening by inserting a sterile probe into the external (secondary) opening and navigating it through the fistula tract. Alternatively, injection of appropriate sterile fluids (saline or hydrogen peroxide) into the external (secondary) opening of the fistula tract and identifying the site of emergence at the internal (primary) opening can also be used to assist in identifying the location of the internal opening.

**NOTE:** Failure to locate the internal (primary) opening may lead to persistence of the fistula. If the internal (primary) opening cannot be reliably identified, an alternative method of treatment should be considered.
3. Using a syringe and short catheter, gently clean the fistula tract as thoroughly as possible by flushing the tract with hydrogen peroxide or sterile saline, or debride using a brush or similar instrument. Do not enlarge the tract.
- 4a. If a seton is not in place, insert a fistula probe, seton, or suitable instrument through the fistula tract, entering through the external (secondary) opening and exiting via the internal (primary) opening. Attach the suture connected to the plug to the instrument at the internal (primary) opening end.
- 4b. If a seton is in place, cut the seton and attach the suture to the seton at the internal (primary) opening end.
5. Draw the tie attached to the plug into the internal (primary) opening and through the fistula tract until slight resistance is felt and the plug securely blocks the internal (primary) opening. Do not use the entire plug unless the tract is the full length of the plug.


**IMPORTANT:** The internal opening is the high-pressure zone of the fistula, as well as the site of ingress of fecal debris. The thicker end of the plug must therefore be securely snugged into the internal (primary) opening to prevent ingress of fecal debris. In addition, the higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract by simple mechanical force.
6. When the plug is properly positioned, trim away and discard any remaining portion of the plug that is not implanted within the fistula tract.

**IMPORTANT:** The plug should be trimmed at the level of the bowel wall at the internal opening in order to minimize contact with bowel contents.
7. Suture the internal (primary) end of the plug in place with suitable resorbable suture. Refer to the illustrations. Suture the internal end of the plug securely to the adjacent tissue, obtaining adequate bites of bowel wall and fistula plug to prevent leakage of bowel contents into the fistula tract and to anchor the fistula plug to prevent migration through the tract.

**NOTE:** There should be no part of the plug visible at the internal (primary) opening.

**NOTE:** If the internal (primary) opening is dimpled or recessed, consider limited mobilization of the mucosal edges prior to suture placement to ensure adequate coverage of plug.

**NOTE:** Inadequate fixation of the plug to the deep tissue layers may result in early recurrence of the fistula.

8. The excess external plug should be trimmed flush with the skin without fixation. The external (secondary) opening may be enlarged to facilitate drainage.

**NOTE:** Complete obstruction of the external (secondary) opening may result in accumulation of fluid, infection, or abscess.

9. Place a sterile dressing over the implant site.
10. Discard any unused portions of the device following standard technique for disposal of medical waste.

#### POST-OPERATIVE CARE

To provide the best environment for tissue integration into the plug, patient activity should be minimized. Provide patients with a list of post-procedure care recommendations. The following patient guidelines should be considered.

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least six weeks.
2. Patients should avoid any heavy lifting over 10lbs (5kg) for at least six weeks.
3. Patients should abstain from sexual intercourse and other forms of vaginal/rectal insertion (e.g. tampons) for at least six weeks.
4. Patients should shower standing up and bathe the area with water to soothe and keep it clean.
5. Patients should observe a liquid diet for 48 hours after surgery.
6. Patients should begin a high fiber diet after 48 hours.
7. Patients should use a stool softener for at least two weeks after surgery.
8. Patients should use the bathroom at least once a day.
9. Patients should take over-the-counter pain medication as needed.
10. Patients should avoid the use of topical steroids.

#### REFERENCES

1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44 (4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection, 1999.* Centers for Disease Control and Prevention (CDC).

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### REFERENCE

- Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.
- Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44 (4):75-80.
- Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection, 1999.* Centers for Disease Control and Prevention (CDC).

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## LITTERATUR

- Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.
- Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44(4):75-80.
- Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection, 1999.* Centers for Disease Control and Prevention (CDC).

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1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.

2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44 (4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

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## REFERENCIAS

1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44 (4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

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1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons.* Clin Infect Dis 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations.* J Med Dent Sci 1997;Dec44 (4):75-80.

3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

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## BIBLIOGRAFIA

1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons*. Clin Infect Dis 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations*. J Med Dent Sci 1997;Dec44 (4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

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#### LITERATUUR

- Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons*. Clin Infect Dis 1997;Apr24(4):609-19.
- Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations*. J Med Dent Sci 1997;Dec44 (4):75-80.
- Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

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## REFERENSER

- Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons*. Clin Infect Dis 1997;Apr24(4):609-19.
- Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations*. J Med Dent Sci 1997;Dec44 (4):75-80.
- Minogram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).

**Keep dry**

**Manufacturer**

**Use-by date**

**Attention, see instructions for use**

**Do not re-use**

**Temperature limit**

**Sterilized using ethylene oxide**

**Date of manufacture**

**MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.  
Phone: 812 339-2235  
Toll Free: 800 457-4500  
Toll Free Fax: 800 554-8335

**EC REPRESENTATIVE**

COOK IRELAND  
O'Halloran Road  
National Technological Park  
Limerick, IRELAND  
Phone: +353 61 334440

**AUSTRALIAN SPONSOR**

WILLIAM A. COOK  
AUSTRALIA PTY. LTD.  
Brisbane Technology Park  
95 Brandl Street  
Brisbane, QLD 4113 Australia  
Phone: +61 7 38 41 11 88