
**COOK**  
MEDICAL

CE 0123

- EN**  
**3** **Biodesign®**  
**Peyronie's Repair Graft**  
Instructions for Use
- CS**  
**4** **Biodesign®**  
**Graft pro reparaci Peyronieovy nemoci**  
Pokyny k použití
- DA**  
**5** **Biodesign®**  
**implantat til reparation af Peyronies sygdom**  
Brugsanvisning
- DE**  
**6** **Biodesign®**  
**Peyronie-Reparaturimplantat**  
Gebrauchsanweisung
- EL**  
**8** **Biodesign®**  
**Μόσχευμα αποκατάστασης της νόσου του Peyronie**  
Οδηγίες χρήσης
- ES**  
**9** **Biodesign®**  
**Injerto para reparación de la enfermedad de Peyronie**  
Instrucciones de uso
- FR**  
**10** **Biodesign®**  
**Grefe pour la réparation de la maladie de la Peyronie**  
Mode d'emploi
- HR**  
**12** **Transplantat Biodesign®**  
**za liječenje Peyroniejeve bolesti**  
Upute za upotrebu
- HU**  
**13** **Biodesign®**  
**graft a Peyronie-kór kezelésére**  
Használati utasítás
- IT**  
**14** **Biodesign®**  
**Innesto per la riparazione per malattia di Peyronie**  
Istruzioni per l'uso
- NL**  
**16** **Biodesign®**  
**Peyronie-reparatieprothese**  
Gebruiksaanwijzing
- NO**  
**17** **Biodesign®**  
**Peyronies reparasjonsimplantat**  
Bruksanvisning
- PL**  
**18** **Biodesign®**  
**wszczep korekcyjny w chorobie Peyroniego**  
Instrukcja użycia
- PT-BR**  
**20** **Biodesign®**  
**Enxerto para reparo da doença de Peyronie**  
Instruções de uso
- SV**  
**21** **Biodesign®**  
**Peyronies reparationsgraft**  
Bruksanvisning
- TR**  
**22** **Biodesign®**  
**Peyronie Onarım Grefti**  
Kullanma Talimatı

\*FP0046-01J\*

## BIODESIGN® PEYRONIE'S REPAIR GRAFT

### INTENDED USE

The Cook® Biodesign® Peyronie's Repair Graft is intended for implantation to reinforce soft tissue where weakness exists in the urological anatomy, including but not limited to repair of tunica albuginea defects and reinforcement in the repair of Peyronie's disease. The graft is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**PEYRONIE'S REPAIR GRAFT** This symbol means the following: Peyronie's Repair Graft

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

### PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- Graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard graft if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that the graft is rehydrated prior to suturing.
- Graft performance has not been evaluated with suture spacing greater than 2 mm.
- Ensure that all layers of the graft are secured when suturing.

### POTENTIAL COMPLICATIONS

The following complications are possible when surgical graft materials are used in the penis. If any of these conditions occur and cannot be managed medically, careful graft removal should be considered:

- Infection
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

**NOTE:** Graft failure resulting in a lack of soft tissue reinforcement may lead to reoccurrence of the Peyronies plaque. If this occurs, additional surgical intervention to remove the plaque may be necessary.

**NOTE:** *De novo* post-operative erectile dysfunction is a known complication associated with reconstructive procedures in the penis. If post-operative erectile dysfunction occurs, graft removal may exacerbate erectile changes. Therefore, graft removal is not recommended. Alternative medical and surgical treatments should be considered.

### STORAGE

This graft should be stored in a clean, dry location at room temperature.

### STERILIZATION

This graft has been sterilized with ethylene oxide.

### INSTRUCTIONS FOR USE

#### Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 50 mL of room temperature sterile saline or sterile lactated Ringer's solution for each graft sheet.

**NOTE:** Handle graft sheets using aseptic technique, minimizing contact with latex gloves.

1. Using aseptic technique, remove the graft inner pouch from its outer bag, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
3. Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
4. Add enough rehydration fluid to the dish to immerse the graft sheet.
5. Allow graft sheet to rehydrate for at least ten minutes.
6. Prepare the graft site using standard surgical techniques.

- Using aseptic technique, trim the graft sheet to fit the site. User feedback suggests that the graft should be sized approximately 30% larger than the site.

**NOTE:** An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to suturing it into place. See step 5.

- Using aseptic technique, transfer the graft sheet to the graft site and suture into place using resorbable suture, while avoiding excess tension.

**NOTE:** Surgical experience indicates that suturing graft sheets with close tissue approximation produces better outcomes.

- Complete the standard surgical procedure.
- Discard any unused portions of the graft sheet following standard technique for disposal of medical waste.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.