COOK
MEDICAL

EN
3

**Biodesign®
Plastic Surgery Matrix**

Instructions for Use

ES
4

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
BIODESIGN® PLASTIC SURGERY MATRIX

INTENDED USE

The Biodesign® Plastic Surgery Matrix is for implantation to reinforce soft tissue where weakness exists in patients requiring soft tissue repair or reinforcement in plastic and reconstructive surgery. The matrix is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

PLASTIC SURGERY MATRIX This symbol means the following:

Plastic Surgery Matrix

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This matrix is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- Matrix is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard matrix if mishandling has caused possible damage or contamination, or if the matrix is past its expiration date.
- Ensure that the matrix is rehydrated prior to cutting or suturing.
- Matrix performance has not been evaluated with suture spacing greater than 2mm.
- Ensure that all layers of the matrix are secured when suturing.
- Suturing more than one matrix together may decrease matrix performance.
- No studies have been done to evaluate the reproductive impact with the clinical use of the matrix.
- The matrix may not have sufficient strength to support stresses encountered in some ventral hernias or large-area, body-wall repairs.
- Patients undergoing radiation therapy may not experience normal wound healing.
- No studies have been done to examine the effect of concomitant radiation with matrix placement.
- Failure to securely suture the matrix to healthy, well-vascularized tissue may result in lack of incorporation of matrix.

NOTE: When implanting the matrix in patients with thin skin, transient, post-operative redness may be seen.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials. If any of these complications occur and cannot be resolved, consider the removal of the matrix:

- Infection
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation.)
- Allergic reaction
- Seroma formation

STORAGE

This matrix should be stored in a clean, dry location at room temperature.

STERILIZATION

This matrix has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle matrix using aseptic technique, minimizing contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
 - Sterile forceps
 - Rehydration fluid: at least 50 mL of room temperature sterile saline or sterile lactated Ringer's solution for each sheet
1. Using aseptic technique, remove the inner pouch containing the matrix from the outer pouch, and place it in the sterile field.

2. Open the inner pouch carefully, and aseptically remove the matrix using sterile forceps.
3. Place the matrix into a sterile dish in the sterile field.
4. Add at least 50 mL of the hydration fluid to the dish.
5. Allow the matrix to re-hydrate until desired handling characteristics are achieved. Rehydration longer than 1 minute is not required and not recommended to exceed 4 minutes.

IMPORTANT: Minimize manipulation of the matrix during rehydration to avoid delamination.

6. Prepare the graft site using standard surgical techniques.
7. Using aseptic technique, trim the matrix to fit the site, providing a small allowance for overlap.

NOTE: An alternative method is to cut the matrix to size prior to rehydration. If this method is selected, be sure to re-hydrate the matrix prior to suturing it into place. See step 5.

8. Using aseptic technique, transfer the matrix to the graft site. Suture into place, ensuring that all layers are captured and avoiding excess tension.

IMPORTANT: Securely suture the matrix to maximize contact with healthy, well-vascularized tissue in order to encourage cell in-growth.

9. Complete the standard surgical procedure.

IMPORTANT: The liberal use of drains is recommended until output is less than 15cc in 24 hours.

10. Discard any unused portions of the matrix according to institutional guidelines for medical waste.

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© COOK BIOTECH 2017
West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.